

EN Titanium Forceps

DE Titan-Pinzetten

ES Pinzas de titanio

Page 115 to 116



ATRAUMA



TITANIUM




OI 10-2401-00
Vena
Closing pressure
110 +/- 10 gr



OI 10-2402-00
Vena
Closing pressure
110 +/- 10 gr



OI 10-2403-00
Vena
Closing pressure
110 +/- 10 gr



OI 10-2404-00
Vena
Closing pressure
110 +/- 10 gr

	Arteria	Arteria	Arteria	Arteria
Closing Pressure	220 +/- 10 gr	220 +/- 10 gr	220 +/- 10 gr	220 +/- 10 gr
Length	30 mm	35 mm	40 mm	50 mm
Jaw Length	9 mm	13 mm	14 mm	23 mm
Weight	1 gr	1 gr	1 gr	1.5 gr



OI 10-2405-00
Vena
Closing pressure
110 +/- 10 gr



OI 10-2406-00
Vena
Closing pressure
110 +/- 10 gr



OI 10-2407-00
Vena
Closing pressure
110 +/- 10 gr



OI 10-2408-00
Vena
Closing pressure
110 +/- 10 gr

	Arteria	Arteria	Arteria	Arteria
Closing Pressure	220 +/- 10 gr	220 +/- 10 gr	220 +/- 10 gr	220 +/- 10 gr
Length	30 mm	35 mm	40 mm	50 mm
Jaw Length	9 mm	13 mm	14 mm	23 mm
Weight	1 gr	1 gr	1 gr	1.5 gr



OI 10-2409-00
Vena
Closing pressure
110 +/- 10 gr



OI 10-2410-00
Vena
Closing pressure
110 +/- 10 gr



OI 10-2411-00
Vena
Closing pressure
110 +/- 10 gr



OI 10-2412-00
Vena
Closing pressure
110 +/- 10 gr

	Arteria	Arteria	Arteria	Arteria
Closing Pressure	220 +/- 10 gr	220 +/- 10 gr	220 +/- 10 gr	220 +/- 10 gr
Length	30 mm	35 mm	40 mm	50 mm
Jaw Length	9 mm	13 mm	14 mm	23 mm
Weight	1 gr	1 gr	1 gr	1.5 gr



ATRAUMA



OI 10-2413-16
16 cm/6 1/4"
1.5 mm



OI 10-2413-20
20.5 cm/8 1/8"
1.5 mm



OI 10-2413-24
24.5 cm/9 5/8"
1.5 mm



OI 10-2414-16
16 cm/6 1/4"
1.8 mm



OI 10-2414-20
20.5 cm/8 1/8"
1.8 mm



OI 10-2414-24
24.5 cm/9 5/8"
1.8 mm

DEBAKEY

